



PREPARE



RESPOND



RECOVER

WESTCOT 15 TIL' 50 PACKAGE

WC15TIL50P

WESTCOT SHELTER IN PLACE W/ IV POLE - WCSIPI (X48): *EXTRUDED ALUMINUM WITH ZINC PLATED STEEL COMPONENTS *FRICTION WELDED TUBULAR ENDS *2" CONVOLUTED FOAM ANTIMICROBIAL SOIL RESISTANT MATTRESS *5 POSITION HEAD ELEVATION WITH LOCKING CORD *ANTI-TILTS LEG DESIGN *BERRY AMENDMENT (10 U.S.C. 2533A) COMPLIANT *ALL MATERIALS EXCEED FED STD 1633 FOR FLAMMABILITY *SAFETY & INSTRUCTIONAL LABELING *FOLDING DIMENSIONS 34" X 42" X 7" *24 LBS *SLEEPING DIMENSIONS 80" X 32" X 18" *WORKING CAPACITY 450 LBS *24 TRIAGE RED *24 TRIAGE YELLOW

WESTCOT DISPOSABLE BULK PACK FITTED SHEET AND BLANKET SET - WCDLPAKFTBT (48) *PREPACKAGED FITTED CHUCK AND BLANKET INSERTED DIRECTLY INTO THE FOLDED WESTCOT. ITEMS ARE DISPOSABLE, INTENDED FOR SINGLE USE

WESTCOT RAPID DEPLOYMENT BAG SHELTER IN PLACE - WCRDSIP2 (X24) *HEAVY DUTY 18 OZ. CUSTOM SEWN VINYL BAG WITH REINFORCED VINYL SKID PLATE *8 SAFETY GRAB HANDLES *WATERPROOF ZIPPER *2 WESTCOTS PER BAG *12 TRIAGE RED *12 TRIAGE YELLOW

WESTCOT ALL TERRAIN CART FULLY LOADED - WCRAT8FL (X6) *HEAVY DUTY WELDED & POWDER COATED STEEL FRAME *4 SHELVES HOLD 2-3 WESTCOTS *5" ALL TERRAIN CASTERS WITH X-AXIS ROTATION & HEAVY DUTY LOCKS *RIVITED IN PLACE HDPE SHELF BASES WITH STEEL SUPPORT BARS *42"X36"X74" @ 500 LBS. (LOADED) *NO ASSEMBLY REQUIRED *OPTIONAL FORK POCKETS

WESTCOT TRIPLE LAYER VINYL DUST COVER - WCVDCWCRAT8 (X6) *ANTIMICROBIAL SOIL RESISTANT DUST COVER FOR WCRAT8FL *SEWN AND SECURABLE ACCESS DOOR *SECURING GROMMETS AND ZIP TIES INCLUDED

WESTCOT TRIAGE CHAIR - WCTR (X20) - *EXTRUDED ALUMINUM WITH ZINC PLATED STEEL COMPONENTS *2 X DECON READY ARM REST *REMOVABLE / REPLACEABLE CHAIR CHAIR COVER *20 TRIAGE GREEN

WESTCOT ALL TERRAIN TRIAGE CHAIR CART - WCATTCC (X2) *HEAVY DUTY WELDED & POWDER COATED STEEL FRAME *HOLDS 10 WCTC *5" ALL TERRAIN CASTERS WITH X-AXIS ROTATION & HEAVY DUTY LOCKS *RIVITED IN PLACE HDPE SHELF BASE WITH STEEL SUPPORT BARS *27X34X58" @ 244 LBS. (LOADED) *NO ASSEMBLY REQUIRED *OPTIONAL FORK POCKETS

WESTCOT TRIPLE LAYER VINYL DUST COVER - WCVDCWCATTCC *ANTIMICROBIAL SOIL RESISTANT DUST COVER FOR WCRAT8FL *SEWN AND SECURABLE ACCESS DOOR *SECURING GROMMETS AND ZIP TIES INCLUDED

XBOMB STANDBY SURGE TRAILER - XBSST - PLUG AND PLAY THIS LEGEND ALL ALUMINUM FLAT NOSE CLIMATE CONTROLLED TRAILER. UNIT COMES WITH 110V POWERED FANS TO CIRCULATE AIR CONTINUOUSLY WHILE IN THE STORAGE AREA.

INCLUDED *48 X WCSIPI *24 X WCRDBSIP2 (12 RED + 12 YELLOW) *48 X WCDLPAKFTBT *6 X WCRAT8FL *6 X WCVDCWCRAT8 *20 X WCTR *2 X WCATTCC *1 X WCVDCWCATTCC *XBSST

"REDUCE RESPONDER INJURY, IMPROVE DEPLOYMENT EFFICIENCY, AND PROTECT YOUR WESTCOT INVESTMENT."