

WESTCOT FUNCTIONAL NEEDS COT

WCFNC

SPECIFICATIONS:

Functional Needs Cot (WCFNC)

Dimensions: 42" x 34" x 8" (Stored)

Dimensions: 80" x 32" x 18" (Deployed)

Dimensions: 42" x 34" x 13.75" (2 Per Carton)

Dimensions: 40" x 48" x 78" @ 500 Lbs (12 Per Pallet)

Weight: 32 Lbs.

Working Capacity: 450 Lbs.

IMS Settle Load Testing: 15,000 Cycles (x 500 Lbs.)

Friction Welded (Spun) Tubular Ends

Heavy Duty Leg Locks

Reinforced Vinyl Coated Decking

Wide Patient Restraint Belt

Antimicrobial Soil Resistant Triple Layer Vinyl Mattress

Antimicrobial Soil Resistant Decking Material

All Materials Exceed FED STD 1633 for Flammability

ADA Compliant

Safety and Instructional Labeling

5 Position Head Elevator with Locking Cord

2 Position Trendelenburg Foot Elevator

2" Convoluted Foam Mattress

Customizable Sleeping Surface

Wheelchair Transferrable Height

Spring and Deck Design

Combined Spring Stretch Resistance: 5,472 Lbs.

Lightweight Rust Proof Aluminum Frame

Stores Securely for Rapid Deployment

Optional Stow on Board Safety Rails

Optional Stow On Board IV Pole

Optional Non-Slip Foot Wraps

Optional Disposable Linen Packs

Optional Rapid Deployment Bag

Optional Push Kit: WCATFNC8P

Application: Functional Needs, Mass Care, Surge Capacity,
MCI, Special Needs, & Medical Shelter

Warranty: Workmanship Warranty & Performance

Guarantee




WESTCOT™