

WESTCOT ALL TERRAIN FUNCTIONAL NEEDS COT 5 PERSON

WCATFNC5P



WESTCOT FUNCTIONAL NEEDS COT SAFETY RAILS + IV POLE - WCFNCSRI (X5):

*EXTRUDED ALUMINUM WITH ZINC PLATED STEEL COMPONENTS *FRICTION WELDED TUBULAR ENDS *2" CONVOLUTED FOAM ANTIMICROBIAL SOIL RESISTANT MATTRESS *5 POSITION HEAD ELEVATION WITH LOCKING CORD *TWO POSITION FOOT ELEVATION *ANTI-TILT *PATIENT SAFETY BELT * STOW ON BOARD SAFETY RAILS & IV POLE *ADA COMPLIANT *BERRY AMENDMENT (10 U.S.C. 2533A) COMPLIANT *ALL MATERIALS EXCEED FED STD 1633 FOR FLAMMABILITY *SAFETY & INSTRUCTIONAL LABELING *FOLDING DIMENSIONS 34" X 42" X 8" *26 LBS *SLEEPING DIMENSIONS 80" X 32" X 18" *WORKING CAPACITY 450 LBS

WESTCOT RAPID DEPLOYMENT BAG FUNCTIONAL NEEDS COTS – WCRDBFNC1

(X5) *HEAVY DUTY 18 OZ. CUSTOM SEWN VINYL BAG WITH REINFORCED VINYL SKID PLATE *5 SAFETY GRAB HANDLES *WATERPROOF ZIPPER *1 WESTCOTS PER BAG

WESTCOT DISPOSABLE LINEN SETS - WCDLPAK (X5) *FITTED SHEET *FLAT SHEET

*BLANKET, *PILLOW *PILLOWCASE *SETS SEALED IN BAG *SET STORED INSIDE OF EACH WESTCOT

WESTCOT ALL TERRAIN CART FULLY LOADED – WCRAT5FL (X1) *HEAVY DUTY

WELDED & POWDER COATED STEEL FRAME *4 SHELVES HOLD 2-3 WESTCOTS *5" ALL TERRAIN CASTERS WITH X-AXIS ROTATION & HEAVY-DUTY LOCKS *RIVITED IN PLACE HDPE SHELF BASES WITH STEEL SUPPORT BARS *42"X36"X74" @ 500 LBS. (LOADED) *NO ASSEMBLY REQUIRED *OPTIONAL FORK POCKETS

WESTCOT TRIPLE LAYER VINYL DUST COVER – WCVDCWCRAT5 (X1)

*ANTIMICROBIAL SOIL RESISTANT DUST COVER FOR WCRAT8FL *SEWN AND SECURABLE ACCESS DOOR *SECURING GROMMETS AND ZIP TIES INCLUDED

INCLUDED *5 X WCFNCSRI *5 X WCRDBFNC1 *5 X WCDLPAK *1 X WCRAT5FL *1 X WCVDCWCRAT5



“REDUCE RESPONDER INJURY, IMPROVE DEPLOYMENT EFFECINCY, AND PROTECT YOUR WESTCOT INVESTMENT.”