

WESTCOT ALL TERRAIN FUNCTIONAL NEEDS COT 8 PERSON

WCATFNC8P



WESTCOT™

WESTCOT FUNCTIONAL NEEDS COT SAFETY RAILS + IV POLE - WCFNCSRI (X8):

*EXTRUDED ALUMINUM WITH ZINC PLATED STEEL COMPONENTS *FRICTION WELDED TUBULAR ENDS *2" CONVOLUTED FOAM ANTIMICROBIAL SOIL RESISTANT MATTRESS *5 POSITION HEAD ELEVATION WITH LOCKING CORD *TWO POSITION FOOT ELEVATION *ANTI-TILT *PATIENT SAFETY BELT *STOW ON BOARD SAFETY RAILS & IV POLE *ADA COMPLIANT *BERRY AMENDMENT (10 U.S.C. 2533A) COMPLIANT *ALL MATERIALS EXCEED FED STD 1633 FOR FLAMMABILITY *SAFETY & INSTRUCTIONAL LABELING *FOLDING DIMENSIONS 34" X 42" X 8" *32 LBS *SLEEPING DIMENSIONS 80" X 32" X 18" *WORKING CAPACITY 450 LBS

WESTCOT RAPID DEPLOYMENT BAG FUNCTIONAL NEEDS COTS - WCRDBFNC2

(X4) *HEAVY DUTY 18 OZ. CUSTOM SEWN VINYL BAG WITH REINFORCED VINYL SKID PLATE *8 SAFETY GRAB HANDLES *WATERPROOF ZIPPER *2 WESTCOTS PER BAG

WESTCOT DISPOSABLE LINEN SETS - WCDLPAK (X8)

*FITTED SHEET *FLAT SHEET *BLANKET, *PILLOW *PILLOWCASE *SETS SEALED IN BAG *SET STORED INSIDE OF EACH WESTCOT

WESTCOT ALL TERRAIN CART FULLY LOADED - WCRAT8FL (X1)

*HEAVY DUTY WELDED & POWDER COATED STEEL FRAME *4 SHELVES HOLD 2-3 WESTCOTS *5" ALL TERRAIN CASTERS WITH X-AXIS ROTATION & HEAVY-DUTY LOCKS *RIVETED IN PLACE HDPE SHELF BASES WITH STEEL SUPPORT BARS *42"X36"X74" @ 500 LBS. (LOADED) *NO ASSEMBLY REQUIRED *OPTIONAL FORK POCKETS

WESTCOT TRIPLE LAYER VINYL DUST COVER - WCVDCWCRAT8 (X1)

*ANTIMICROBIAL SOIL RESISTANT DUST COVER FOR WCRAT8FL *SEWN AND SECURABLE ACCESS DOOR *SECURING GROMMETS AND ZIP TIES INCLUDED

INCLUDED *8 X WCFNCSRI *4 X WCRDBFNC2 *8 X WCDLPAK *1 X WCRAT8FL *1 X WCVDCWCRAT8



WESTCOT™

"REDUCE RESPONDER INJURY, IMPROVE DEPLOYMENT EFFECINCY, AND PROTECT YOUR WESTCOT INVESTMENT."