

WESTCOT ACTIVE PATIENT CARE WCAPC

SPECIFICATIONS:

Dimensions: 42" x 34" x 9" (Stored)
Dimensions: 80" x 32" x 28" (Deployed)
Dimensions: 42" x 35" x 10" (1 Per Carton)
Dimensions: 40" x 48" x 78" @ 500 Lbs (10 / Per Pallet)
Weight: 36 Lbs.
Working Capacity: 450 Lbs.
IMS Settle Load Testing: 15,000 Cycles (x 500 Lbs.)
Friction Welded (Spun) Tubular Ends
Heavy Duty Leg Locks
Reinforced Vinyl Coated Decking
Wide Patient Restraint Belt
Antimicrobial Soil Resistant Triple Layer Vinyl Mattress
Antimicrobial Soil Resistant Decking Material
All Materials Exceed FED STD 1633 for Flammability
Safety and Instructional Labeling
5 Position Head Elevator with Locking Cord
2 Position Trendelenburg Foot Elevator
2" Convoluted Foam Mattress
Customizable Sleeping Surface
Wheelchair Transferrable Height
Spring and Deck Design
Combined Spring Stretch Resistance: 5,472 Lbs.
Lightweight Rust Proof Aluminum Frame
Stores Securely for Rapid Deployment
Optional Stow on Board Safety Rails
Optional Stow On Board IV Pole
Optional Non-Slip Foot Wraps
Optional Disposable Linen Packs
Optional Rapid Deployment Bag
Optional Push Kit: WCATAPC6P
Application: Field Hospital, Surge Capacity, MCI, Special
Needs, & Medical Shelter
Warranty: Workmanship Warranty & Performance
Guarantee




WESTCOT™