

DonorEase™ Mobile

For Portability and Efficiency in Operations

- Lightweight Aircraft Aluminum Frame
- Replaceable Mattress Deck
- Integrated Cushioned Seats
- Fully Adjustable Head & Foot Sections from Lounge to Trendelenberg
- Interchangeable Molded Arm Rest
- Dual Integrated Trip Scale Mounts
- Multi-position Pillow
- Fastened Sub-Assembly Construction for Ease of In-House Maintenance
- Weight Capacity - 400 lbs.

Options and Accessories

- Trip Scale Rod
- Head, Foot and Arm Rest Covers

Dimensions

- 10" D x 28" W x 40" T - Folded
- 81" L x 28" W x 24" T - Open
- 33 lbs. - Bed and Accessories



Contact us about
Design Features,
Options and Purchasing

Shaun Taylor
855.505.9267 or 704.466.3874
sales@integritymedicalsolutions.com

IMS4BEDS.COM



Made in USA