



IMS4BEDS.COM
855.505.9267



Functional Needs Cot (WCFNC)

Dimensions

42 x 32 x 8" (Stored)
80 x 32 x 18" (Deployed)
30 Lb. Curb Weight
450 Lb. Working Capacity
42 x 34 x 13.75 (Carton / Stores 2)

Safety

Exceeds IMS Settle Load Testing at 15,000 Cycles (x 450 Lbs.)
Friction Welded (Spun) Tubular Ends
Heavy Duty Leg Locks
Reinforced Vinyl Coated Decking
Wide Patient Restraint Belt
Antimicrobial Soil Resistant Triple Layer Vinyl Mattress
Antimicrobial Soil Resistant Decking Material
Easy to Decontaminate
All Materials Exceed FED STD 1633 for Flammability
ADA Compliant
Safety and Instructional Labeling

Comfort & Client Benefits

5 Position Head Elevator with Locking Cord
2 Position Trendelenburg Foot Elevator
2" Convoluted Foam Mattress
Customizable Sleeping Surface
Wheelchair Transferrable Height
Spring and Deck Design
Combined Spring Stretch Resistance of 5,472 Lbs.

Convenience

Lightweight Rust Proof Aluminum Frame
Stores Securely for Rapid Deployment
Optional Stow on Board Safety Rails
Optional Stow On Board IV Pole

Application

Functional Needs Care, Mass Care, Surge Capacity, Hospital Surge Capacity Programs, Special Needs, & Medical Shelter

Warranty

Workmanship Warranty & Performance Guarantee



W
E
S
T
C
A
R
T

O
P
T
I
O
N

Tube Enterprises, Inc. dba Integrity Medical Solutions, Inc. TID #: 20-1590871 / 20-4604067
855.505.9267 SALES@INTEGRITYMEDICALSOLUTIONS.COM IMS4BEDS.COM
DUNS # 144934929 CAGE CODE # 8A0Q5