



IMS4BEDS.COM
855.505.9267



Hot Cot 24 (WCHC24)

Dimensions

42 x 24 x 8" (Stored)
80 x 24 x 18" (Deployed)
30 Lb. Curb Weight
450 Lb. Working Capacity

Safety

Exceeds IMS Settle Load Testing at 15,000 Cycles (x 450 Lbs.)
Friction Welded (Spun) Tubular Ends
Heavy Duty Leg Locks
Reinforced Vinyl Coated Decking
Wide Patient Restraint Belt
Antimicrobial Soil Resistant Triple Layer Vinyl Mattress
Antimicrobial Soil Resistant Decking Material
Easy to Decontaminate
All Materials Exceed FED STD 1633 for Flammability
ADA Compliant
Safety and Instructional Labeling

Comfort & Client Benefits

5 Position Head Elevator with Locking Cord
2 Position Trendelenburg Foot Elevator
2" Convoluted Foam Mattress
Customizable Sleeping Surface
Wheelchair Transferrable Height
Spring and Deck Design
Combined Spring Stretch Resistance of 5,472 Lbs.

Convenience

Lightweight Rust Proof Aluminum Frame
Stores Securely in EMS Vehicles and in small compartments
Optional Stow on Board Safety Rails
Optional Stow On Board IV Pole

Application

EMS Vehicle Rapid Response, Surge Capacity, Special Needs, Mass Care, Functional Needs Care, & Medical

Warranty

IMS products have a standard One Year Warranty on Workmanship and Product Performance. IMS' Managing Partners stand behind the Westcot Brand and Guarantee Complete Performance Satisfaction.

