



IMS4BEDS.COM
855.505.9267



Bariatric Patient Care W/Head Elevator (WCBPCHE)

Dimensions

42 x 38 x 7" (Stored)
82 x 38 x 22" (Deployed)
42 Lb. Curb Weight
700 Lb. Working Capacity

Safety

Exceeds IMS Settle Load Testing at 15,000 Cycles (x 700 Lbs.)
Friction Welded (Spun) Tubular Ends
Heavy Duty Leg Locks
Reinforced Vinyl Coated Decking
Wide Patient Restraint Belt
Antimicrobial Soil Resistant Triple Layer Vinyl Mattress
Antimicrobial Soil Resistant Decking Material
Easy to Decontaminate
All Materials Exceed FED STD 1633 for Flammability
ADA Compliant
Safety and Instructional Labeling

Comfort & Client Benefits

2" Convuluted Foam Mattress
Heavy Duty Head Elevator
Customizable Sleeping Surface
Easy Entry Height of 22"
Spring and Deck Design
Combined Spring Stretch Resistance of 9,576 Lbs.

Convenience

Lightweight Rust Proof Aluminum Frame
Stores Securely for Rapid Deployment
Elevate Client's Upper Body
Optional Safety Rails
Optional Stow On Board IV Pole

Application

Above Average Size and/or Weight Shelter Clients, Pregnant Woman
Comfort, Special Needs, Functional Needs, Mass Care, Surge Capacity,
Hospital Surge Capacity Programs, & Field Medical Hospital

Warranty

IMS products have a standard One Year Warranty on Workmanship and
Product Performance. IMS' Managing Partners stand behind the
Westcot Brand and Guarantee Complete Performance Satisfaction.

