



IMS4BEDS.COM
855.505.9267



Bariatric Patient Care (WCBPC)

Dimensions

42 x 38 x 6" (Stored)
82 x 38 x 22" (Deployed)
36 Lb. Curb Weight
700 Lb. Working Capacity

Safety

Exceeds IMS Settle Load Testing at 15,000 Cycles (x 700 Lbs.)
Friction Welded (Spun) Tubular Ends
Heavy Duty Leg Locks
Reinforced Vinyl Coated Decking
Wide Patient Restraint Belt
Antimicrobial Soil Resistant Triple Layer Vinyl Mattress
Antimicrobial Soil Resistant Decking Material
Easy to Decontaminate
All Materials Exceed FED STD 1633 for Flammability
Safety and Instructional Labeling
Optional Anti-Slip Feet

Comfort & Client Benefits

2" Convoluted Foam Mattress
Customizable Sleeping Surface
Easy Entry Height of 22"
Spring and Deck Design
Combined Spring Stretch Resistance of 9,576 Lbs.

Convenience

Lightweight Rust Proof Aluminum Frame
Stores Securely for Rapid Deployment
Optional Stow On Board IV Pole

Application

Above Average Size and/or Weight Shelter Clients, Pregnant Woman
Comfort, Functional Needs Care, Mass Care, Surge Capacity, Hospital
Surge Capacity Programs, & Field Medical Hospital

Warranty

IMS products have a standard One Year Warranty on Workmanship
and Product Performance. IMS' Managing Partners stand behind the
Westcot Brand and Guarantee Complete Performance Satisfaction.

