



IMS4BEDS.COM
855.505.9267



Patient Transport Unit (WCPTU)

Dimensions

- 42 x 28 x 10" (Stored)
- 80 x 28 x 28" (Deployed)
- 44 Lb. Curb Weight
- 350 Lb. Working Capacity

Safety

- Exceeds IMS Settle Load Testing at 15,000 Cycles (x 450 Lbs.)
- Friction Welded (Spun) Tubular Ends
- Heavy Duty Leg Locks
- Reinforced Vinyl Coated Decking
- Wide Patient Restraint Belt
- Antimicrobial Soil Resistant Triple Layer Vinyl Mattress
- Antimicrobial Soil Resistant Decking Material
- Easy to Decontaminate
- All Materials Exceed FED STD 1633 for Flammability
- Safety and Instructional Labeling

Comfort & Client Benefits

- 5 Position Head Elevator with Locking Cord
- 2 Position Trendelenburg Foot Elevator
- 2" Convuluted Foam Mattress
- Customizable Sleeping Surface
- Spring and Deck Design
- Combined Spring Stretch Resistance of 5,472 Lbs.

Convenience


- X-Axis 360 Rotation (2 Swivel/Locking, 2 Sraight, 2 Swivel)
- Lightweight Rust Proof Aluminum Frame
- Stores Securely for Rapid Deployment
- Optional Stow on Board Safety Rails
- Optional Stow On Board IV Pole

Application

Indoor Use Only: Mass Care, Surge Capacity, Hospital Surge Capacity Programs, Special Needs Transport, & Field Medical Hospital

Warranty

IMS products have a standard One Year Warranty on Workmanship and Product Performance. IMS' Managing Partners stand behind the Westcot Brand and Guarantee Complete Performance Satisfaction.



W
E
S
T
C
O
R
T

O
P
T
I
O
N