

General Use Cot (GUC)

Dimensions

42 x 32 x 6" (Stored)
80 x 32 x 22" (Deployed)
22 Lb. Curb Weight
450 Lb. Working Capacity

Safety

Exceeds IMS Settle Load Testing at 15,000 Cycles (x 450 Lbs.)
Friction Welded (Spun) Tubular Ends
Heavy Duty Leg Locks
Reinforced Vinyl Coated Decking
Antimicrobial Soil Resistant Triple Layer Vinyl Mattress
Antimicrobial Soil Resistant Decking Material
Easy to Decontaminate
All Materials Exceed FED STD 1633 for Flammability
Safety and Instructional Labeling

Comfort & Client Benefits

2" Convoluted Foam Mattress
Customizable Sleeping Surface
Wheelchair Transferable Height – 18"
Spring and Deck Design
Combined Spring Stretch Resistance of 5,472 Lbs.

Convenience

Lightweight Rust Proof Aluminum Frame
Stores Securely for Rapid Deployment

Application

Mass Care, Surge Capacity Programs, Functional Populations,
Homeless Sheltering Operations, Budget Friendly

Warranty

IMS products have a standard One Year Warranty on Workmanship
and Product Performance. IMS' Managing Partners stand behind the
Westcot Brand and Guarantee Complete Performance Satisfaction.



W
E
S
T
C
O
T



O
P
T
I
O
N