

Westcot All Terrain Functional Needs Package 8 Person WCATFNCP8



The Westcot All-Terrain Functional Needs Cot 8 Person Package is a Rapid Deployment System designed to assist Response Agencies with the deployment of essential Functional Needs Cots and equipment to nearly any environment. The Westcart AT8 heavy-duty deployment cart outfitted with four foam filled pneumatic all-terrain casters which allows for fluid movement and zero deflation concerns. Two of the all-terrain casters afford high strength locking systems which increase operational safety. There are eight Westcot Functional Needs Cots (WCFNC) paired inside of four Westcot Rapid Deployment Bags (WCRDBW2). The WCFNCs are designed to provide equal care for Special Needs Survivors of Natural Disasters by providing functional head and leg elevation as well as stow on board IV poles, Safety Rails, and a set of Disposable Linen packed inside of each WCFNC. The WCRDBW2 were designed to reduce Responder injury, increase deployment efficiency, and protect your Westcot investment. The whole package is covered in a Triple Layer Vinyl which is antimicrobial, soil resistant, and puncture resistant.

WCAT8 Dimensions / Total Loaded Weight: 39 x 43 x 74 / 466 Lbs.

WCFNCsRI Dimensions / Weight – 80 x 32 x 18 (Capacity 450 Lbs.) / 40 x 34 x 8 (33 Lbs. Each)

WCRDBW2 – Fits Two WCFNCsRI – enables proper lifting and rapid deployment

WCDLP – 1 Fitted Sheet, 1 Flat Sheet, 1 Blanket, 1 Pillow Case, 1 Pillow packed inside of each cot