



IMS4BEDS.COM

855.505.9267



Bariatric Special Needs 5 Person System (WCBSN5)

Includes:

- Qty 5 WCBSN
- Qty 5 WCBSN Protective Covers
- Qty 15 Disposable Linen Packs
- Qty 1 WCRBSN5
- Qty 1 Vinyl Dust Cover

Dimensions:

- 42 x 38 x 7" (Stored)
- 82 x 38 x 22" (Deployed)
- 42 Lb. Curb Weight

Safety

- 700 Lb. Working Capacity
- Exceeds IMS Settle Load Testing at 15,000 Cycles (x 700 Lbs.)
- Friction Welded (Spun) Tubular Ends
- Heavy Duty Leg Locks
- Reinforced Vinyl Coated Decking
- Wide Patient Restraint Belt
- Antimicrobial Soil Resistant Triple Layer Vinyl Mattress
- Antimicrobial Soil Resistant Decking Material - Easy Decon.!
- All Materials Exceed FED STD 1633 for Flammability
- Safety and Instructional Labeling
- ADA Compliant

Comfort & Client Benefits

- 2" Convoluted Foam Mattress
- Heavy Duty Head Elevator
- Customizable Sleeping Surface
- Easy Entry Height of 22"
- Spring and Deck Design
- Combined Spring Stretch Resistance of 9,576 Lbs.

Convenience

- Lightweight Rust Proof Aluminum Frame
- Stores Securely for Rapid Deployment
- Elevate Client's Upper Body
- Optional Safety Rails
- Optional Stow On Board IV Pole

Application

- Above Average Size and/or Weight Shelter Clients, Pregnant Woman
- Comfort, Special Needs, Functional Needs, Mass Care, Surge Capacity, Hospital Surge Capacity Programs, & Field Medical Hospital

Warranty

IMS products have a standard One Year Warranty on Workmanship and a Product Performance Guarantee.

