

Westcot Active Patient Care 8 Person WCAPC8P



The Westcot Active Patient Care 8 Person (WCAPC8P) offering is a rapid deployment system for WCAPC which were designed to provide operating height (28") to patients and/or medical professional tasked with emergency response functions. The WCAPC8P was created to improve hospital surge capacity programs or agencies tasked with field medial hospital activation.

There are eight Active Patient Care Westcots (WCAPC) paired inside of four Westcot Rapid Deployment Bags (WCRDBAPC2) stored inside of one Westcart 8a (WC8A.) The stow on board IV poles and Safety rails provide necessary precaution, the Disposable Linen Packs (WCDLP) improves deployment and decontamination efficiencies. The Whole system is covered in a Triple Layer Vinyl Dust Cover (WCVDC) which is soil and puncture resistant. This system reduces Responder injury, increase deployment efficiencies, and protects your Westcot investment.

WC8A: 39 x 43 x 82" at 155 lbs.

WCFNCSRIV: 80 x 32 x 28" (deployed) / 40 x 34 x 9" at 35 lbs. (stored)

WCRDBAPC2: Holds 2 x WCAPCSRIV

WCDLP: 1 Fitted Sheet, 1 Flat Sheet, 1 Blanket, 1 Pillowcase, 1 Pillow, stored inside of a Ziploc bag