

# • DonorEase™ Tilt Features

## Product Features

- Powder Coated Tubular Steel Frame w/Ergonomic Design
- Stainless Steel Tilt Control Bar, Gas Spring, w/Quick Release Lever
- Medical Grade Upholstery - Easy to Clean, Antimicrobial, Stain Resistant, Abrasion Resistant, Multiple Colors available
- Multi-Position Upholstered Pillow
- Solid Rubber Molded Arm System w/Locking T-Handle and R/L Stanchions.
- Optional Plug & Play Comes w/Two 120V Receptacles and Two USB Ports
- Integrated “Donor Cubby” for Storage of Donor’s Personal Effects
- Adjustable Leg Levelers
- Mattresses and Arm Rests are Interchangeable between Power, Tilt, and Comfort Chairs



# • Features Identified

