

• DonorEase™ Comfort Option

Product Features

- One Piece Powder Coated Tubular Steel Frame w/Ergonomic Design
- Medical Grade Upholstery – Easy to Clean, Antimicrobial, Stain Resistant, Abrasion Resistant, Multiple Colors available
- Multi-Position Upholstered Pillow
- Solid Rubber Molded Arm System w/Locking T-Handle and R/L Stanchions.
- Optional Two 120V Receptacle and Two USB Ports
- Integrated “Donor Cubby” for Storage of Donor’s Personal Effects
- Adjustable Leg Levelers
- Mattresses and Arm Rests are Interchangeable between Power, Tilt, and Comfort Chairs




Comfort Wedge Pillow Solution



The new DonorEase™ EZ-T™ (Trendelenburg) Pillow allows you to elevate your donor's legs quickly and easily for immediate treatment of a reaction. The EZ-T™ can be stationed throughout the center for quick retrieval. We recommend a bright easy to see color, so that the EZ-T™ can quickly be located and sprung into action.

The EZ-T™ is made from the same high quality Naugahyde™ as the mattress and head pillow, cleaning and decontamination are simple and effective.

The pillow elevates from approximately 3.5" to 11". At 20" wide, it supports your donor's legs and rests safely on the DonorEase™ chair.



• Comparative Specifications

	<u>Tilt</u>	<u>Comfort</u>
Total Weight	107 lbs.	107 lbs.
Entry Height	24"	24"
Length	68"	68"
Length Tilted	72"	N/A
Width, Overall	29"	29"
Useable Mattress Width	24"	24"
Height Foot End	31"	31"
Height Head End	47"	47"
Weight Capacity	600 lbs.	600 lbs.



VS

