

• DonorEase™ Comfort Option

Product Features

- One Piece Powder Coated Tubular Steel Frame w/Ergonomic Design
- Medical Grade Upholstery – Easy to Clean, Antimicrobial, Stain Resistant, Abrasion Resistant, Multiple Colors available
- Multi-Position Upholstered Pillow
- Solid Rubber Molded Arm System w/Locking T-Handle and R/L Stanchions.
- Optional Two 120V Receptacle and Two USB Ports
- Integrated “Donor Cubby” for Storage of Donor’s Personal Effects
- Adjustable Leg Levelers
- Mattresses and Arm Rests are Interchangeable between Power, Tilt, and Comfort Chairs

